


PULSE LAVAGE SYSTEM



Disposable Pulse Lavage System for Cleaning procedures in modern Cementing Techniques

Thorough and effective cleaning. Used for hip/knee joint implantations, trauma and wound debridement.

Optimum cementing fixation and penetration results as well as reduced chance of blood lamination.

The use of the Pulse Lavage System will prevent microembolization of the marrow contents and minimise circulatory changes. It is also proved that the Pulse Lavage System reduces the risk of fat embolism.

Technical Specifications:

Weight: 210 g

Nozzle Pattern: Elliptical Orifice

Flow (l/min): 1000/200

Powered by: 8 x 1,5 V Cell

System is supplied with a high capacity fan spray debridement tip and narrow intra modulary tip!

Art. No. PL 1000

infection
guaranteed